## QUICK GUIDE

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v 1.0 2021

# QUICK GUIDE

ELECTRONIC PRECISION BALANCES (LCD, GRAPHIC AND TOUCHSCREEN DISPLAY)


This is just a guide for quick start using the instrument. The complete **USER MANUAL** describing in depth correct use and all functions, parameters, and features of the instrument must be downloaded from the link reported below.

Use of the instrument different from the one reported in the user manual does not ensure anymore the safety of the product.

## GET READY TO WORK

- **Remove** the balance and the accessories from the package and check if the instrument has not suffered any damage due to transportation. Do not throw away the original package.
  - **In models with internal calibration** remove the brass screw, if present, in the left lower part of the balance (see “balance parts” chapter of this guide)
  - **Do not install** the balance in places with air flows, heavy thermal changes and vibrations.
  - **Place** the support pan and weighing pan on the balance (**Note**: 0,01mg balances don't have support pan)
  - **Level the balances** using the levelling feet and looking at the level bubble until is centered
  - **Insert** the supplied power supply unit jack into the connector on the back of the instrument then connect the power adapter to a power outlet near the instrument; do not use cables/extensions that do not comply with applicable regulations. Check that the power shown on the device plate label corresponds to the one in use in the country where you are installing the item.
  - **Balance Warm-up time:** wait for 8-12 hours from switch-on for Semimicro (0,01mg) and Analytical balances (0,1mg) and 30 minutes for other models. Then calibrate the balance using an appropriate weight following the instructions of the paragraph “calibration” of this guide.
- For a correct use of the instrument never disconnect the balance from the mains. If you want to**

turn it off, use the **ON / OFF** button to put it in **Stand-by** (in this way it is not necessary to wait each time the warm-up time).



- **Service** must be performed by specialized personnel and the spare parts to use must be genuine. Always contact the dealer from whom the instrument has been purchased.

#### **Important precautions for 0,01mg and 0,1mg models**

- Place the balance in a secluded location: not put it near doors to prevent air drafts. Avoid places with many persons. Protect the balance from an air conditioner or ventilation fans.
- Keep the **room temperature** as constant as possible, at a value between 15 and 30°C.
- Keep the **humidity** as constant as possible, it should be between 40% and 65%.
- Place the balance away from heat sources, eg heaters, light bulbs (use tubular lamps), windows (the warmth of the sun's rays can filter out the window and may affect the weighing result).
- The table where the balance is placed must be **stable** (eg. Laboratory bench or marble / granite countertop): must not sag during lean of the instrument.
- Must be as **anti-magnetic** and **anti-static** as possible. Must be **reserved** to the balance.
- Place the balance as close as possible to the table legs because there are less **vibrations** rather than in the middle of the table.
- Load the sample to weigh always in the **center of the plate** to avoid possible errors.
- Open the glass doors of the **cabinet** as less as possible and use the **tweezers** to load/unload the samples to be weighed.
- Pay attention to possible **static charges** that may be generated due to containers with materials of low electrical conductivity or due to dry air (less than 40% moisture). Electrostatic charges can alter the results of the weighings. It is recommended to use a ionizer (optional) to eliminate static charges present on samples or accumulated on balance's parts.  
*Because of static charges the results of weighings are always different, the balance is not repeatable.*
- Pay attention to **dynamic push**: a big **temperature** difference between the sample to be weighed and the weighing chamber creates air drafts along the sample. A colder object appears heavier while a warmer object lighter, this effect is reduced when it is reached the thermal equilibrium between sample and weighing chamber.  
*With the dynamic push you will get results that move in one direction depending on whether the material is colder or warmer.*
- Pay attention to substances that can **evaporate** (alcohol) or **absorb humidity** (silicon gel). Because of these types of materials weight may vary constantly in one direction.
- Be careful to **magnetic** materials: magnetic objects will attract each other, the force that results are wrongly interpreted as a load.  
*With magnetic materials the weighing results are hardly repeatable, the indication remains stable but weighing result provides different results.*

# ASSEMBLY


Model series: from SM to G


Balances semimicro with resolution **0,01mg** have support pan and weighing pan in one single piece, as shown in the photo.

Model series: H, I, L, M




Model series: HC, HCG. TS-HC



## BALANCE PARTS

|    |                                                                                                                                                                                        |
|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | Weighing pan                                                                                                                                                                           |
| 2  | 9 poles (pin) female connector for RS232 interface for connection to printer/PC (not for M models)                                                                                     |
| 3  | Levelling feet                                                                                                                                                                         |
| 4  | connector for power supply                                                                                                                                                             |
| 5  | Screws for balance closing. In models with 4 feet it is necessary first to remove the two rear feet (one fixed and one moving) to access to the closing screws                         |
| 6  | Hook for under-balance weighing. Hook not available in internal calibration models. <b>Note:</b> please remove the cap to access the hook                                              |
| 7  | Protection screw for autocalibration motor (only for models with internal calibration). <b>ATTENTION: remove this screw by hand as soon as the balance is removed from its package</b> |
| 8  | Reset button (Press with a sharp object inside the hole to reset the instrument when necessary)                                                                                        |
| 9  | USB port (only in Touchscreen models)                                                                                                                                                  |
| 10 | Switch on/off key (only in Touchscreen models)                                                                                                                                         |




### Model series: from SM to G

### Model series: H, I, L, M

### Model series: HC, HCG, TS-HC

# KEYBOARD and DISPLAY SYMBOLS

## LCD and GRAPHIC DISPLAY models

| FIGURE                                                                              | BUTTON           | FUNCTION                                                                                                                                            |
|-------------------------------------------------------------------------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MENU             | SHORT KEY PRESSING:<br>call-up application menu of the balance<br><br>LONGTIME PRESSING<br>UNTIL SOUND STOPS:<br>enter/exit of user menu of balance |
|                                                                                     | Arrow down       | menu scroll forward arrow                                                                                                                           |
|    | ON/OFF switch    | Standby button (OFF) or switch on (ON) . Exit button (ESC)                                                                                          |
|    | TARE and Zeroing | TARING, ZEROING                                                                                                                                     |
|                                                                                     | Arrow up         | left scroll arrow                                                                                                                                   |
|  | PRINT            | SEND data to printer/interface                                                                                                                      |
|                                                                                     | Arrow ENTER      | CONFIRM selection                                                                                                                                   |
|  | CALIBRATION      | balance CALIBRATION                                                                                                                                 |
|                                                                                     | Arrow up         | menu scroll backward arrow                                                                                                                          |

| SYMBOL DISPLAYED                                                                    | MEANING                 |
|-------------------------------------------------------------------------------------|-------------------------|
| *                                                                                   | Stability indicator     |
| O                                                                                   | Zero indicator          |
| %                                                                                   | Percentage weighing     |
| PC                                                                                  | Piece counting          |
|  | Battery level indicator |
| ▼                                                                                   | Data input mode         |
| H                                                                                   | Higher threshold        |
| L                                                                                   | Lower threshold         |
| DS                                                                                  | Density measure         |
| ct, ozt, lb, GN, dwt, Kg, mg                                                        | Measuring unit          |

## TOUCH SCREEN models

## 1 - Indication area for weighing values

|                                                                                   |                                                                                                                        |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>g</b>                                                                          | Measuring unit indication, by tapping the symbol, you can access the measurement unit menu and select the desired unit |
| *                                                                                 | Stability indicator                                                                                                    |
| <b>T Man</b>                                                                      | Indication of weight displayed. <b>Net</b> net value, <b>Gro</b> gross value                                           |
| <b>T =</b>                                                                        | Tare value indicator                                                                                                   |
|  | Analog indicator bar of magnitude being measured.                                                                      |
| Max 6200g d=0.01g                                                                 | Indication of the metrological data of the balance                                                                     |
| <b>T = 12.81 g</b>                                                                | Information about tare value                                                                                           |

## 2 - Bar displaying date/time, user selection, access to functions, settings, device shutdown.

|                                                                                    |                                                                                                                                        |
|------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
|   | Tap this area to access the available <b>functions</b> menu.                                                                           |
|   | Tap this area to access the menus for weighing <b>settings</b> , display preferences, peripheral units setup, and users customization. |
|  | Tap this area to access the <b>switching</b> on or off menu for the device standby mode.                                               |
| Default                                                                            | Tap this area to access the <b>users</b> menu to quickly select the one desired.                                                       |

How to enter the **calibration** menu of the balance:

|   |                    |              |
|---|--------------------|--------------|
| < | Calibration Menu   | ×            |
|   | Calibration        |              |
|   | Calibration data   |              |
|   | Calibration mode   | Ext. default |
|   | Calibration weight | 2000 g       |

Press the Settings menu button on the main window and select the "**Calibration menu**"

In the calibration menu, you can calibrate, display and print the data for the last calibration and set the calibration mode. Before performing the calibration procedure, check that the plate is empty. By tapping the calibration key, the calibration procedure corresponding to the set calibration mode will start. Load (in case of external calibration) the required calibration weight and wait for weight acquisition. If the operation is successful, the weighing screen will display with the calibrated weight value. Remove the weight from the plate.

## 3 - By tapping the keys displayed in this area, the relative function will be executed. Keys available in simple weighing mode:


|       |       |     |     |                                                                                                                                                                                                                                                 |
|-------|-------|-----|-----|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Print | Gro   | -T- | -0- | <b>Print:</b> Data printing activation key.<br><b>-0-:</b> Zero activation key.<br><b>-T-:</b> Tare operation activation key.<br><b>Gro/Net:</b> Net or gross weight display mode key.<br><b>T Man:</b> Key for entering a tare value manually. |
| Print | T Man | -T- | -0- |                                                                                                                                                                                                                                                 |

## SWITCH-ON

In “**STANDBY**” mode press the **ON/OFF** button (in verified balance models, after switch-on of the balance the display visualizes the message “WARM UP” and a countdown). Place goods to be weighed on the pan, wait the stability symbol ✱ (asterisk) to appear and then read the weight value. Press again **ON/OFF** button to return to “**STANDBY**” mode.

## LEVELLING

Level the balance with level bubble, adjusting the Front Feet until the air bubble is placed in the centre of the indicator. To lift the balance → rotate front feet clockwise. To Lower the balance → rotate feet anticlockwise. Once level bubble is properly placed, lock the feet screwing the disk foot-locking (not in all models)

## TARING

Place the container on the pan and press **O/T** button. “**O-t**” is displayed. As soon as the stability is reached the value of zero “**0.000**” will be displayed (weight of container is internally stored). Place the goods to be weighed in the container, the net weight will now be displayed.

## CALIBRATION of the balance

For accurate measurements, the balance must be adjusted to the local environment, using the calibration function. Select the desired calibration mode in the Calibration menu of the balance.

**NOTE: if there is any interference or instability during calibration process, an error message will be visualized.**

**External calibration** (external calibration is not available in verified balance models)

Press the **CAL** button when no object is on the weighing pan. Wait until the value of required calibration weight appears flashing on the display, then put the required calibration weight carefully in the center of the pan. Display will stop flashing. After calibration is performed, remove the calibration weight. Balance automatically returns to weighing mode.

If you want to calibrate with a different (must be a greater multiple) weight from the calibration one, long press the **CAL** button to start the LOAD function. Calibration then proceeds the same as described above.

**Internal calibration** (ONLY for models with internal calibration)

Press the **CAL** button when the weighing pan is unloaded. You will see the message “**CAL**” displayed and then the balance calibration will be automatically performed. At the end of the calibration the balance returns automatically to standard weighing mode

**Automatic calibration** (ONLY for models with internal calibration)

In this mode the balance performs an automatic calibration when the temperature variation exceeds a (factory set) value or at time intervals (factory set), the first of the two that occurs first. The autocalibration is performed through the internal mass and only when the weighing pan is unloaded and only when no other weighing operations are being performed. At the end of the calibration the balance returns automatically to standard weighing mode.

If, due to vibrations or strong air-flows, the calibration is not terminated, the message “**CAL but**” is displayed. Press again the **CAL** button to retry.

If the problem does not disappear, select the “external calibration” mode and perform a calibration, then contact the service center.

In this mode it is also possible to perform a calibration with the internal mass by pressing the **CAL** key manually at any time, being sure that weighing pan is unloaded.

# ACCESS TO BALANCE FUNCIONS and SETTINGS MENU

## LCD and GRAPHIC DISPLAY MODELS

Press tue **MENU** button to enter the balance functions menu (Piece-counting, Percentage waighing,...), then browse the menu through **MENU** and **CAL** buttons and confirm choices using the **PRINT** button.

Long press the **MENU** button to enter the balance configuration menu (Measuring units, Language, Date setup,...) of the balance, then browse the menu through **MENU** and **CAL** buttons and confirm choices using the **PRINT** button.

## TOUCHSCREEN DISPLAY MODELS

Tap this area to access the available functions menu.

Tap this area to access the menus for weighing settings, display preferences, peripheral units setup, and users customization.

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**OPTIKA® S.r.l.**