# OPTIKA®

B A L A N C E S  
I T A L Y

## QUICK GUIDE

Model
SM Series
SMG Series
A Series
AG Series

v 1.0 2019

# OPTIKA®

B A L A N C E S  
I T A L Y

EC-Konformitätserklärung  
EC-Déclaration de conformité  
EC-Dichiarazione di conformità

EC-Declarção de conformidade  
EC-Declaration of conformity  
EC-Declaración de conformidad

OPTIKA S.r.l., Via Rigla 30 Ponteranica 24010 (BG) ITALY

GB	Declaration of conformity	We hereby declare that the product to which this declaration refers complies with the following standards.
D	Konformitätserklärung	Wir erklären hiermit, dass das Produkt, auf das sich diese Erklärung bezieht, mit den nachstehenden Normen übereinstimmt.
F	Déclaration de conformité	Nous déclarons avec cela responsabilité que le produit, auquel se rapporte la présente déclaration, est conforme aux normes citées ci-après.
I	Dichiarazione di conformità	Dichiariamo con ciò che il prodotto al quale la presente dichiarazione si riferisce è conforme alle norme di seguito citate.
E	Declarção de conformidade	Manifestamos en la presente que el product al que se refiere esta declaración está de acuerdo con las normas siguientes..
P	Declaración de conformidad	Declaramos por meio da presente que o produto no qual se refere esta declaração, corresponde às normas seguintes

## Electronic Balance: Series SM,SMG, A, AG

European Directive	Harmonized Standard
2014/30/EU (EMC)	EN 61326-1:2013
2014/35/EU (LVD)	EN 61010-1:2010 EN 61010-2-010:2014
2011/65/EU (RoHS)	EN 50581:2012

This declaration of conformity is issued under the sole responsibility of the manufacturer.

# QUICK GUIDE

## ELECTRONIC PRECISION BALANCES

Series: SM, SMG, A, AG

### ATTENTION:

This is a guide to start quickly working and do NOT substitute the complete user manual; we please you to read completely the user manual for a correct installation and a right use of your new electronic balance.

Usage of the instrument different from the one reported in the user manual does not ensure anymore the safety of the product.

## 1 Positioning the instrument

### 1.1 Preface

Balances with resolutions 0.1mg and 0.01mg are instruments of high sensitivity and precision. To achieve reliable and accurate results it is fundamental to install the instrument in a suitable environment to meet the requirements necessary to ensure its proper functioning.

### 1.2 Place recommended to install the instrument

Choice of the correct placement of the instrument is fundamental to ensure optimal and precise operations.

**For optimum weighing it is necessary to respect the following criteria:**

#### - TYPE OF ENVIRONMENT

- Place the balance in a corner of the room to minimize vibrations.
- Place the balance in a secluded location: not put it near doors to prevent air drafts.
- Avoid places with many persons: every person who is close generates a draft of air.
- Protect the balance from an air conditioner or ventilation fans and other electronic devices with fans (eg. Computer or other laboratory instruments).
- Keep the room temperature as constant as possible, at a value between 15 and 30 ° C. You must use a conditioner (but leave the speed of the fan to a minimum to avoid excessive air drafts).
- Keep the humidity where balance is used as constant as possible, it should be between 40% and 65% of moisture.
- Place the balance away from heat sources, eg heaters, light bulbs (use tubular lamps), windows (the warmth of the sun's rays can filter out the window and may affect the weighing result).

## - TYPE TABLE

- The table must be **stable**: must not sag during lean of the instrument (eg. Laboratory bench or marble / granite counter top).
- Must be as **anti-magnetic** and **anti-static** as possible.
- Must be **reserved** to the balance.
- Place the balance as close as possible to the table legs because there are less **vibrations** rather than in the middle of the table.

## 2 Recommendations for a correct use of the instrument

### 2.1 Preface

To obtain accurate and repeatable weighing pay attention to the following.

### 2.2 First power on

- **THE FIRST TIME THAT YOU CONNECT THE BALANCE TO THE MAINS, BEFORE USE, WAIT AT LEAST 12 HOURS TO ENSURE THE WARM-UP OF THE INSTRUMENT.**
- For a correct use of the instrument never disconnect the balance from the mains. If you want to turn it off, use the ON / OFF button to put it in **Stand-by** (in this way it is not necessary to wait each time the warm up time).

### 2.3 Using the Balance

- **Level the balance.** It is fundamental always to level properly the instrument: regularly check that the air bubble is at the center of the level. Help with **block-feet** to always ensure the correct leveling of the balance.
- Load the sample to weigh always in the **center of the plate** to avoid possible errors.
- Open the glass doors of the **cabinet** as less as possible and use the **tweezers** to load/unload the samples to be weighed.
- Pay attention to possible **Static charges** that may be generated due to containers with materials of low electrical conductivity or due to dry air with less than 40% moisture. Electrostatic charges can alter the results of the weighings. It is recommended to use the **lonizer Mod. Ion-A15** to eliminate static charges present on any substances to be weighed or that have accumulated on the instrument.
  - Because of static charges the results of weighings are always different, the balance is not repeatable.
- Pay attention to **dynamic push**: a big **temperature** difference between the sample to be weighed and the weighing chamber creates air drafts along the sample. A colder object

appears heavier while a warmer object lighter, this effect is reduced when it is reached the thermal equilibrium between sample and weighing chamber.

- With the dynamic push you will get results that **move in one direction** depending on whether the material is colder or warmer.

- Pay attention to substances that can **Evaporate** (alcohol) or **Absorb humidity** (silicon gel). Because of these types of materials weight may vary **constantly in one direction**.
- Be careful to **Magnetic** materials: magnetic objects will attract each other, the force that results are wrongly interpreted as a load.


- With magnetic materials the weighing results are hardly repeatable, the indication remains stable but weighing result provides different results.

## 2.4 Care and cleaning of the balance

- Before **cleaning**, remove all removable parts (eg. Plate, underplate).
- To clean the weighing chamber, the plate, the glasses and the other parts please use an **antistatic** liquid (do not use cloths that can leave fragments).


### 3 Balance Overview

#### 3.1 Front Balance view

- 1 Draft windshield
- 2 Electronics box
- 3 Function keys
- 4 Antiventilation ring
- 5 Display LCD
- 6 Adjustable front foot

#### 3.2 Rear Balance view



- 7 Level bubble
- 8 Connector 9 poles (pin) female for interface RS232 for printer PC
- 9 Rear fixed foot
- 10 Label with balance model and serial number
- 11 Power socket


# 4 Keyboard and display

Series: **SM, A**


**SMG, AG**


Standby (OFF/ON) or escape ESC.


TARE or zero button..

Selection CONFIRM or SEND data to printer.

Balance setup MENU button.

Balance CALIBRATION button.

\*

Stability indicator

O

Zero indicator

%

Percentage weighing

PC

Piece counting

Battery charge indicator

Insert data mode

H

Upper threshold

L

Lower threshold

DS

Density measure mode

ct, Measure unit

ozt, lb, GN, dwt, Kg, mg


## 5 Operating

### 5.1 Components location


- As a first step place **Pan Support (1)** on balance cone
- Lay **Weighing Pan (2)** on pan support and check the pan is placed properly.
- Put **Antiventilation ring (3)**

### 5.2 Connection balance and electronic box

- Insert plug connection cable into jack of electronic box.
- Fix the connector screwing the two side screws.

### 5.3 Connection between power supply and electronics box


- Connection to network is made through an external power supply.
- Insert power supply connector in DC plug placed on electronics box.
- Use **ONLY** power supply provided with the balance.

## 5.4 Balance levelling

- **Level** the balance with level bubble, adjusting it ( **2** ) until the air Bubble is placed at the centre of the indicator ( **1** ).

- **Adjust** the balance according to the level bubble:

**Lift** the balance → rotate front feet clockwise

**Lower** the balance → rotate feet anticlockwise

- Once level bubble is properly placed using adjusting feet ( **2** ), lock them screwing the disk Foot-locking ( **3** )

## 6 Switch on

After connecting the balance to power supply, an autodiagnosis of electronic circuits is automatically effected, ending with stand by indication

**Balance Warm-up:** Wait 8 hours from switch on for warm up

**It is suggested to never disconnect the balance from power socket and use ON/OFF key to put the instrument in standby mode until end of work.**

From “**STAND BY**” mode: to bring the balance back to working conditions, press **ON/OFF** key.

It is recommended not to drop heavy objects on balance pan, to avoid damage of the instrument.


Electronic balance effects mass measurings using gravity (g). Differences in geographical areas e in altitude change gravity acceleration (g).

Therefore, to get precise measurings, balance has to be adjusted to environmental conditions. This adjusment is accomplished through calibration function.

It is needed to calibrate the balance every time it is moved to another place.

## 7 WEIGHING


In “**STAND BY**” mode press the **ON/OFF** button. Place goods to be weighed on the pan, wait the stability symbol ✱ (asterisk) to appear and then read the weight value

Press again **ON/OFF** button to return to “**STAND BY**” mode.


## 8 TARING

Place the container on the pan and press **O/T** button. “**O-t**” is displayed. As soon as the stability is reached the value of zero “**0.000**” will be displayed (weight of container is internally stored). Place the goods to be weighed in the container, the net weight will now be displayed

## 9 EXTERNAL CALIBRATION and LOAD FUNCTION

Thanks to **LOAD** function of the balance, it is also possible to effect calibration of the balance using a weight with a value greater than the factory setting value:


2

While the **LOAD** message is flashing, load on the pan a weight with value equal or greater to the factory setting calibration weight value; the balance recognizes as valid a weight equal or greater to factory setting calibration weight IF its value is an integer value considering the most significant digit. *For example:* if factory setting value of calibration weight is 200g, then it is possible to calibrate the balance also with weights which value starts from 200g, 300g, 400g up to the upper limit of the balance weighing range.

3

As the calibration is effected with success, the display stops flashing and the value of the calibrated weight is visualized with the current measuring unit. Remove the calibration weight used from the pan of the balance. The balance is ready for weighing.

**NOTE:** if there is any interference during calibration process, an error message will be visualized.

## 10 INTERNAL CALIBRATION


### 10.1 Autocalibration (AUT-CAL)

The balance makes an auto-calibration when the temperature variation exceeds the factory value and/or at factory set time intervals. The autocalibration is performed through the internal mass and only when the weighing plate is unloaded and only when no other weighing operations are being performed.

In this mode it is also possible to perform a calibration with the internal mass by pressing the **CAL**-key manually at any time, being sure that weighing plate is unloaded.


1

Press the **CAL** button when the weighing plate is unloaded. You will see the message “**CAL**” displayed and then the balance calibration will be automatically performed

2

At the end of the calibration the balance returns to standard weighing mode

IF, due to vibrations or strong air-flows, the calibration is not terminated, the message “**CAL bUt**” is displayed. Press again the **CAL** button to retry.

If the problem does not disappear, select the “external calibration” mode and contact the service center.

# 11 OVERVIEW OF BALANCE PARAMETERS MENU


SEQUENCE	1	2	3	4	5
ACTION	<b>MENU BUTTON</b> LONG TIME PRESSING	<b>MENU BUTTON</b>	<b>CAL BUTTON</b>	<b>ENTER BUTTON</b>	<b>MENU BUTTON</b> LONG TIME PRESSING
RESULT	you enter in balance parameters set menu	scroll backward the menu	scroll forward the menu	confirm the choice selected	you exit from balance parameters set menu

<i>Parameters</i>	<i>Description</i>
UNITS	select measuring units
PC-PRINT	select print mode
BAUD RT	set velocity of serial communication
AUTO 0	auto-zero setting
FILTER	filters level
STABIL	stability level
BLT	backlight mode
TIME OFF	auto switch off
CONTR	display contrast
CALIB	calibration mode

For a complete description of functions, it is warmly suggested to read the user manual of the balance

## 12 OVERVIEW OF BALANCE PROGRAMS MENU

Press the **MENU** button. Press again **MENU** button to scroll forward to next parameter or press **CAL** button to return to previous parameter; press **PRINT** button to confirm the choice. To exit the menu, press and keep pressed the **MENU** button until the acoustic signal ends.  
For a complete description of balance programs, it is warmly suggested to read the user manual of the balance

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**OPTIKA® S.r.l.**