OPTIKA[®]

B A L A N C E S
I T A L Y

QUICK GUIDE

Model
TG Series

v 1.0 2019

Declaration of Conformity

EC-Konformitätserklärung
EC-Déclaration de conformité
EC-Dichiarazione di conformità

EC-Declaração de conformidade
EC-Declaration of conformity
EC-Declaración de conformidad

OPTIKA S.r.l., Via Rigla 30 Ponteranica 24010 (BG) ITALY

GB	Declaration of conformity	We hereby declare that the product to which this declaration refers complies with the following standards.
D	Konformitätserklärung	Wir erklären hiermit, dass das Produkt, auf das sich diese Erklärung bezieht, mit den nachstehenden Normen übereinstimmt.
F	Déclaration de conformité	Nous déclarons avec cela responsabilité que le produit, auquel se rapporte la présente déclaration, est conforme aux normes citées ci-après.
I	Dichiarazione di conformità	Dichiariamo con ciò che il prodotto al quale la presente dichiarazione si riferisce è conforme alle norme di seguito citate.
E	Declaración de conformidad	Por la presente, manifestamos que el producto al que se refiere dicha declaración, estan fabricados bajo las siguientes normativas
P	Declaração de conformidade	Declaramos por meio da presente que o produto no qual se refere esta declaração, corresponde às normas seguintes

Electronic Balance: Series TG

European Directive	Harmonized Standard
2014/30/EU (EMC)	EN 61326-1:2013
2014/35/EU (LVD)	EN 61010-1:2010 EN 61010-2-010:2014
2011/65/EU (RoHS)	EN 50581:2012

QUICK GUIDE

ELECTRONIC PRECISION MOISTURE BALANCES


Series: TG

ATTENTION:

This is a guide to start quickly working and do NOT substitute the complete user manual; we please you to read completely the user manual for a correct installation and a right use of your new electronic balance.
Usage of the instrument different from the one reported in the user manual does not ensure anymore the safety of the product.

1. INSTALLATION INSTRUCTIONS

- **Remove** the balance, the calibration weight (if included) and the accessories from the package and check if any visible damage of the instrument
- **Do not install** the balance in places with air flows, heavy thermal changes and vibrations. **Place** the support pan and weighing pan on the balance
- **Level the balances** using the levelling feet and looking at the level bubble 
- **Wait for 30 minutes from switch on** and then calibrate the balance (using the appropriate mass if included) following the instructions of the paragraph "calibration"

WARNING:

Please follow carefully these steps for installing and use the new balance before starting your work routine. A way of use of the instrument different from this user manual will not guarantee the instrument's safety anymore. Keep this manual in a safe place.

Please follow these indications to avoid problems and to grant a safe use of moisture analyzer:

- Use the moisture analyzer exclusively for the moisture determination of samples. Every incorrect use of this instrument may create danger for safety of persons and may cause damages at the instrument or at other objects.
- Do not use the instrument in areas where there is danger of explosions; also make the instrument work only respecting the environment conditions as reported in this manual.
- If this material is used in plants or in environment conditions that require strict safety conditions, please respect the instructions of the directives regarding the installation of this material currently existing in your country.
- This instrument must be used only by qualified personell, who knows the features and characheristics of the sample under test.
- Before starting to work with the instrument for the first time, please check that voltage is the same of your electric mains.
- To get the instrument free from tension, remove the power supply cable.
- Place straight the power supply cable so that to avoid contact with the very hot surfaces of the instrument.
- Use only extension cables that are conform to normatives and that are equipped of a protection conductor.

Attention, protection against heat

- Respect the distance and the free space around the instrument to avoid an accumulation of heat in the instrument and and the over heating of the instrument ifself:
 - 20 cm around the instrument
 - 1 m upon the instrument
- Do not place flammable materials upon, below or near the instrument since the heater warms the surrounding area.
- Remove the samples carefully, the heater and the dishes for samples may be still very hot.

Danger for persons or objects when working with particular samples: Fire Explosion

- Flammable or explosive substances
- Substances that contain solvents
- Substances that during drying emit gas or flammable explosive vapors

Package content

- 1 Device for humidity determination.
- 2 VDE power cable
- 3 15 pin M/F cable for connecting the balance/heater
- 4 Tablemat
- 5 Sample-holder plate pan
- 6 Pan with antiventilation cylinder
- 7 #10 sample-holder plates
- 8 CD with instruction manual

ENGLISH


Start-up

The following describes all the steps required to prepare the device for the first start-up.

Open the oven lid and place:

1. Pan with the antiventilation cylinder
2. Place the sample-holder plate extractor
3. Place the star on the weighing cone.

The balance is connected to the heater via a 15-pin M/F cable.


Plug the cable into the two connectors on the rear of the device as shown in the figure.

Plug the VDE cable supplied into the supply connector on the rear of the device.

N.B. Verify that the power supply indicated on the device nameplate matches the one in use in the country where you are installing the unit.

So connect the VDE cable to the outlet in the vicinity of the device. Do not use cables/extensions non-compliant with current regulations.


Wait 30 minutes after the start-up and calibrate the device prior its leveling. For the calibration procedure, refer to Chapter 7.6. Perform calibration of the device each time it is moved to another place.

It is recommended not to drop overweight items on the balance weighing pan to avoid damaging it.

The service must be performed by trained personnel and spare parts used must be original. To do that, please contact the retailer from whom the device was purchased.

Display area functions

The thermobalance is equipped with a color display with tap screen functionality that allows tapping on the screen in the active areas to quickly perform the access to the various menus and functionality of the device.

1

Area showing the weighing values or data related to humidity during and at the end of the drying cycle

g Display of the measurement unit; tapping the symbol, you can access the menu of the measurement units and select the ones desired.

***** Stability value display.

T = Tare value display.

0%-100% Analog bar showing the magnitude being measured.

2

Area showing the heating operation status

60 °C Value display of the heater temperature.

Stb: 60°C Value display of the standby temperature of the heater.

High temperature danger display.

3

Bar to display the time and date, to set the weighing parameters, to select the user, to access to drying functions and to switch off the device.

By tapping this area, you access the menu of the drying mode functions

By tapping this area, you access the weighting menu settings, the viewing preferences, the devices and user customization setup

By tapping this area, you access the switching off menu or enable the standby status of the device

By tapping this area, you access the user menu to make your selection quickly than desired





4

Bar to display the function keys.

By tapping the keys displayed in this area, the relative function will be executed. The following are all available keys that may vary depending on the selected applications

Standby, logout, switching off functions

After inserting the power cord into the wall outlet, the device will automatically switch on taking in the basic weighing screen.

By tapping the area with the symbol  You access the switching off menu.

There are 3 switching off modes:

- Display off.
- Logout.
- Switching off.

Display off

Selecting the display off mode, the device will go in the standby state; the display is turned off, so that only the bottom right command to turn on the system will be enabled.

Tapping the symbol, you will return to the basic weighing screen.

Logout

Selecting the logout mode, you will disable the current user

Tapping **“Start”**, it will activate the user and in the case of password protected user it will be asked you to enter it.

Tapping  you can select another user by which you want to switch on the device.


Switching off

Selecting the switching off command, the device will be turned off. To turn it on, you shall press the key on the left of the device or remove the power plug and insert it again.


Determination of the humidity content with "Quick" method

This section describes how to perform the analysis of the humidity content using the "Simple" method that allows you to perform a quick test by defining only the temperature and the method of drying end.

To access the function, tap the "Simple" icon in the menu "Functions".


The display screen allows you to start a drying cycle by tapping the "Start" button with the temperature parameters and modalities of the drying end displayed to the right of the function.

- **Temperature.**

Enter the desired temperature and tap the "OK" key.

The temperature values that can be set range from a minimum of 35°C to a maximum of 160°C) with step of 1°C.

- **Drying end.**

The end drying parameter is the setting that defines the method by which you want to determine the end of the humidity measuring cycle contained in the substance to be analyzed.

< 1 ... 99 min X

10 min

7	8	9	←
4	5	6	C
1	2	3	OK
+/-	0	.	

Time: setting this method of the end drying, the cycle will end when the time set will be elapsed. Also you can end the cycle manually at any time by tapping the "Stop" key.

Enter the value of time and press the "OK" button to confirm

< Simple X

Start	
Temperature	100°C
Drying end	Auto 0.10% / 30s

Autostop %: setting this method of the end drying, the cycle will end when the variation of humidity loss will be less than 0.1% for a time interval of 30 seconds. Also you can end the cycle manually at any time by tapping the "Stop" key.

Start: tapping "Start", you will enable the procedure for the humidity content analysis.

Load the test pan on the weighing plate and close the oven. You will automatically tare the test plate.

< Pan weighing

 -L- °C
 3.204* g

0%

 100%

Load the pan and close the heater

Abort -0-

Load the substance to analyze and close the oven.

< Sample weighing

 -L- °C
 1.229 g

0%

 100%

Load the sample and close the heater


Abort -T- -0-

N.B. The minimum amount of substance required for the drying cycle based on the balance resolution:

- 0.01g** Min 500mg
- 0.001g** Min 500mg
- 0.0001g** Min 50mg

After closing the heater cover, the cycle will automatically start and the analysis screen will appear.

During drying, they are displayed all the parameters related to the test and through a real-time graph you can see the trend of drying.

For all the details related to the phase of drying and the final result management, refer to **Chapter 11**.

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OPTIKA® S.r.l.