

# OPTIKA quick finder

Resolution	Weighing Range	Family Name	Page	M Verification	Internal Calibration	External Calibration	Windshield
0,01mg	42g...220g	TS - TSI	8-9		■		■
0,01mg	42g...220g	SM	14		■		■
0,01mg	42g...220g	SMG	15		■		■
0,1mg	220g...610g	TS - TSI	10-11		■		■
0,1mg	410g...610g	A	16		■		■
0,1mg	410g...610g	AG	17		■		■
0,1mg	120g...310g	B	18	□	□	■	■
0,1mg	120g...310g	BG	19	□	□	■	■
0,001g	210g...1000g	TS	12-13		■	■	□
0,001g	520g...1200g	C	20	□	□	■	□
0,001g	520g...1200g	CH	21	□	□	■	■
0,001g	520g...1200g	CG	22	□	□	■	□
0,001g	520g...1200g	CGH	23	□	□	■	■
0,001g	160g...420g	E	26		□	■	□
0,001g	160g...420g	EG	28		□	■	□
0,001g	160g...420g	EH	27		□	■	■
0,001g	160g...420g	EGH	29		□	■	■
0,001g	120g...420g	H	34			■	■
0,001ct / 0,001g - 0,01g	330g...4200g	JW	42		▲	■	▲
0,01g	3100g...6200g	TS	12-13		■	■	
0,01g	2200g...8200g	D	24	□	□	■	
0,01g	2200g...8200g	DG	25	□	□	■	
0,01g	1000g...5200g	F	30	□	□	■	
0,01g	1000g...5200g	FG	31	□	□	■	
0,01g	620g...4200g	I	35			■	
0,01g	200g...420g	M	37			■	
0,01g/0,1g	8kg...20kg	TS	12-13			■	
0,1g	4500g...20000g	G	32			■	
0,1g	4500g...20000g	GG	33			■	
0,1g	3200g...6500g	L	36			■	
0,1g	1000g...2200g	M	37			■	
0,01g-0,1g-1g	6200g...32000g	TS-HC - HCG - HC	43-45			■	

## Moisture Balances

Resolution	Weighing Range	Family Name	Page	Thermometer Calibration	Halogen Infrared Lamp
0,1mg	60g	TG	38-39	□	■
0,1mg	60g	T GRAPHIC	40-41	□	■
0,001g	160g	TG	38-39	□	■
0,001g	160g	T GRAPHIC	40-41	□	■
0,01g	60g	TG	38-39	□	■
0,01g	60g	T GRAPHIC	40-41	□	■

( ■ = Standard □ = Optional ▲ = Only For Some Models )

# OPTIKA quick finder

Lcd Display	Graphic Display	Touch Screen display	Integrated Ioniser	Rs232	Units	Glp	Multi Language	Hook	Battery / Battery Pack
		■	▲	■	■	■	■		
■			▲	■	■				
	■		▲	■	■	■	■		
		■	▲	■	■	■	■		
■			▲	■	■				
	■		▲	■	■	■	■		
■			▲	■	■				
	■		▲	■	■	■	■		
		■		■	■	■	■	▲	
■				■	■			▲	
	■			■	■	■	■	■	
				■	■			▲	□
	■			■	■	■	■	▲	□
■				■	■			▲	□
	■			■	■	■	■	▲	□
■				■	■			■	□
				■	■	■	■	■	□
■				■	■			■	□
	■			■	■	■	■	■	□
■				■	■			■	□
		■		■	■	■	■	■	□
■				■	■			■	□
	■			■	■	■	■	■	□
▲	▲	▲		■	■	▲	▲	■	□

Graphic Display	Touch Screen Display	Rs232	Units	Glp	Multi Language	Internal Database, Save results	Password protection
	■	■	■	■	■	■	■
■		■	■	■	■		
	■	■	■	■	■	■	■
■		■	■	■	■		
	■	■	■	■	■	■	■
■		■	■	■	■		